



Analytical Gas Systems Products for LC/MS & Evaporation

Bulletin AGS-LCMS

aerospace
climate control
electromechanical
filtration
fluid & gas handling
hydraulics
pneumatics
process control
sealing & shielding

Parker Hannifin Corporation

The Global Leader in Motion and Control Technologies

We engineer success of our customers around the world, drawing upon nine core motion and control technologies. These technologies enable virtually every machine and process to operate accurately, efficiently and dependably.

As the global leader in motion and control, we partner with our distributors to increase our customers' productivity and profitability by delivering an unmatched breadth of engineered components and value-added services.

We continue to grow with our customers by creating application-focused products and system solutions. A key to our global expansion has been to follow our customers and establish operations, sales and service wherever they are needed. No single competitor matches Parker's global presence.



Corporate Headquarters
in Cleveland, Ohio.

Parker's Motion and Control Technologies

Aerospace	Hydraulics
Climate Control	Pneumatics
Electromechanical	Process Control
Filtration	Sealing & Shielding
Fluid & Gas Handling	

Legal Notifications



FAILURE OR IMPROPER SELECTION OR IMPROPER USE OF THE PRODUCTS AND/OR SYSTEMS DESCRIBED HEREIN OR RELATED ITEMS CAN CAUSE DEATH, PERSONAL INJURY AND PROPERTY DAMAGE.

This document and other information from Parker Hannifin Corporation, its subsidiaries and authorized distributors provide product and/or system options for further investigation by users having technical expertise. It is important that you analyze all aspects of your application and review the information concerning the product or system in the current product catalog. Due to the variety of operating conditions and applications for these products or systems, the user, through its own analysis and testing, is solely responsible for making the final selection of the products and systems and assuring that all performance, safety and warning requirements of the application are met.

The products described herein, including without limitation, product features, specifications, designs, availability and pricing, are subject to change by Parker Hannifin Corporation and its subsidiaries at any time without notice.

Offer of Sale

The items described in this document are hereby offered for sale by Parker Hannifin Corporation, its subsidiaries or its authorized distributors. This offer and its acceptance are governed by the provisions stated in the "Offer of Sale".

© Copyright 1995, 2012, Parker Hannifin Corporation, All Rights Reserved.

Table of Contents

	Products for LC/MS & Evaporation	Page
	NitroVap Gas Generators	4
	Nitrogen Generators for LC/MS Flow Capacities up to 60 lpm	6
	NitroFlow Lab Self Contained LC/MS Membrane Nitrogen Generator	8
	SOURCE LC/MS TriGasGenerator Series Model LCMS-5000NA	9
	SOURCE LC/MS TriGas Generator Series Model LCMS-5001NTNA	11
	SOURCE LC/MS TriGas Generator Series Model LCMS-5001TNA	13
	SOURCE LC/MS Super Flow TriGas Generator Series Model LCMS-SF5000NA	15
	Low and Mid Flow Nitrogen Generators	17
	High Flow Nitrogen Generators	19
	Voltage Appendix	25
	Recommended Gas Generators	26
	Additional Parker Products	28
	Gas Generator Services	29
	Offer of Sale	31

NitroVap Gas Generators

Parker Balston's NitroVap-1LV and NitroVap-2LV Nitrogen Generators can provide clean, ultra-dry nitrogen to sample evaporators. These systems offer high nitrogen output flows in a miniature cabinet. The user can activate the manual SLEEP economy mode to eliminate compressed air consumption when the sample concentrator is not in use.

Nitrogen is produced by utilizing a combination of filtration and membrane separation technologies. A high efficiency prefiltration system pretreats the compressed air to remove all contaminants down to 0.01 micron. Hollow fiber membranes subsequently separate the clean air into a concentrated nitrogen output stream and an oxygen enriched permeate stream, which is vented from the system. The combination of these technologies produces a continuous on demand supply of pure nitrogen.

The NitroVap generators are complete systems with state-of-the-art, highly reliable components engineered for easy installation, operation, and long term performance. The Parker Balston NitroVap-1LV and NitroVap-2LV eliminate all the inconveniences and cost of LN2 dewar and nitrogen cylinder gas supplies and dependence on outside vendors. Uncontrollable price increases, dewar ice and condensation, contract negotiations, long term commitments, and tank rentals are no longer a concern. With a NitroVap generator, you control your gas supply.

Since NitroVap generators incorporate unique membrane separation technology, nitrogen delivery is immediate to the sample concentrator. "Lock-it-and-leave-it" operation of the sample concentrator is maintained without downtime and without "running out of gas" mid blow-down.



NitroVap-1LV and NitroVap-2LV



This Technology Features Advanced HiFluxx Fiber

Features and Benefits

- Ideal for any combination of sample evaporators up to 100 nozzle positions
- Produces clean, dry (to -20°F) dewpoint evaporator grade nitrogen from any standard laboratory compressed air source
- Accelerates evaporation by decreasing the partial vapor pressure above the solvent liquid
- Eliminates inconvenient and dangerous LN2 boil-off dewars and nitrogen gas cylinders from the laboratory
- Recommended and used by many sample concentrator and sample evaporator manufacturers
- Payback period of typically less than one year
- Sleep economy mode
- Silent operation and minimal operator attention required

NitroVap Gas Generators


Principal Specifications

Nitrogen Purity	Up to 90%
Nitrogen Dewpoint	Down to -20°F (-29°C) atmospheric
Maximum Nitrogen Flow Rate	NitroVap-1LV: up to 80 slpm @ 100 psig input up to 140 slpm @ 125 psig input NitroVap-2LV: up to 160 slpm @ 100 psig input up to 287 slpm @ 125 psig input
Electrical Requirements	None
Nitrogen Outlet Pressure	0-15 psig user controlled
Dimensions	10.63"w x 14.1"d x 16.5"h (26.92cm x 35.81cm x 41.91cm)
Inlet Port/Outlet Port	1/4" NPT (female)
Shipping Weight	53 lbs/24 kg

Ordering Information

for assistance, call 800-343-4048, 8 to 5 Eastern Time

Description	Model
NitroVap Nitrogen Generators	NitroVap-1LV and NitroVap-2LV
Maintenance Kit (Includes 1 each filter cartridge, and 1 each membrane cartridge)	MKNITROVAP
Preventive Maintenance Plan	NITROVAP-1LV-PM, NITROVAP-2LV-PM
Extended Support with 24 Month Warranty	NITROVAP-DN2



Use with These and Other Blowdown Evaporators

- TurboVap from Biotage
- N-Evap from Organomation
- RapidVap from LabConco
- Reacti-Vap from Fisher Pierce
- Duo-Vap from Jones Chromatography
- DryVap from Horizon Technology
- Evaporex from Apricot

Self Contained Nitrogen Generators for LC/MS Flow Capacities up to 60 lpm

The Parker Balston® NitroFlow 60 is a self contained generator that produces up to 60 slpm of pure LC/MS grade nitrogen at pressures of up to 110 psig. Nitrogen is produced utilizing a combination of a scroll compressor and nitrogen membrane separation technologies. This combination of technologies yields the highest performing, most reliable and quietest integrated nitrogen generation system available.

The NitroFlow 60 is also available with an integrated membrane dryer for use with instruments that require dry air, including the chip cube interface from Agilent Technologies.

Typical applications include LC/MS, nebulizer gases for APCI and ESI, Jet Stream, I Funnel, ELSD, Turbo Vaps and chemical solvent evaporation.

The unique combination of a rotary scroll compressor and high efficiency membrane ensures that the NitroFlow 60 has many unique advantages over all other existing LC/MS nitrogen generators. Rotary scroll compressors operate at low temperatures, have less moving parts and are significantly quieter than piston compressors used by other Nitrogen Generator manufacturers.



Features and Benefits

- Complete “plug and play” system recommended for all major LC/MS instruments
- Phthalate-free, no organic vapors
- Produces a continuous supply of nitrogen for all LC/MS applications
- Eliminate dangerous nitrogen cylinders from the laboratory
- Nearly silent operation; operates at less than 49 dB(A)

Nitrogen Generators for LC/MS

Flow Capacities up to 60 lpm


Principal Specifications

Model	NitroFlow 60	NitroFlow 60D
Nitrogen	Up to 60 slpm	Up to 60 slpm
Dry Air Flow	N/A	5 slpm
Dry Air Dewpoint	N/A	-40°F (-40°C)
Hydrocarbon Free	Yes	Yes
Phthalate Free	Yes	Yes
Maximum Outlet Pressure	100 psig	100 psig
Atmospheric Dewpoint	-58°F (-50°C)	-58°F (-50°C)
Outlet Port	Female 1/4" NPT	Female 1/4" NPT
Min/Max Ambient Temperature	50°F/95°F (10°C/35°C)	50°F/95°F (10°C/35°C)
Electrical Requirements	195-254 VAC, 60 Hz, 1 Phase, 14A* 230 VAC, 50 Hz, 1 Phase, 13A*	195-254 VAC, 60 Hz, 1 Phase, 14A* 230 VAC, 50 Hz, 1 Phase, 13A*
Dimensions	43"H x 21"W x 34"D (109cm H x 53cm W x 86cm D)	43"H x 21"W x 34"D (109cm H x 53cm W x 86cm D)
Shipping Weight	643 lbs. (292 kg)	643 lbs. (292 kg)

* During operation, 30A at startup

Ordering Information for assistance, call 800-343-4048, 8 to 5 Eastern Time

Description	Model Number
Nitrogen Generator with Integrated Compressor	NitroFlow 60
Nitrogen Generator with Dryer Option and Integrated Compressor	NitroFlow 60D
NitroFlow 60-PM	PM service program
NitroFlow 60-PMPLUS	Plus service program
NitroFlow60-DN2	Depot extended warranty
NitroFlow60-EN2	Express extended warranty



NitroFlow Lab Self Contained LC/MS Membrane Nitrogen Generator

The Parker Balston® LC/MS NitroFlow Lab is a self-contained membrane nitrogen generator that produces LC/MS grade nitrogen with output pressure to 116 psig. Nitrogen is produced by utilizing a combination of compressors, carefully matched with filtration, and membrane separation technology components.

Intake ambient air from the laboratory is filtered using an inlet suction breather filter to remove airborne organic and particulate impurities. This purified air is delivered to a long life low pressure air compressor which provides an air stream to hollow fiber membranes which subsequently separate the clean air into a concentrated nitrogen retentate and

oxygen enriched permeate, which is then cycled through the system. Prior to exiting the system pure nitrogen retentate is delivered to a nitrogen amplification compressor to assure proper pressure, flow and purity to the LC/MS. The Parker Balston LC/MS NitroFlow Lab will deliver a continuous or on demand supply of pure nitrogen making it the smart alternative to cylinders. Superior engineering with carefully matched filtration, membrane separation and compression technologies have resulted in a system with the utmost reliability and longevity. Additional applications include: nebulizer gases, chemical and solvent evaporation, instrument supply and purge, evaporative light scattering equipment and sparging.



NitroFlow Lab



Features and Benefits

- Flow capacity to 30 LPM
- Includes 2 year compressor warranty
- Ideal for all derivatives of ESI and APCi modes
- Includes state-of-the-art, oil-less compressors
- Unlike PSA Hosmer technologies, membrane will not suppress corona needle discharge
- Special sound insulation design ensures quiet operation

Principal Specifications

Model	NitroFlowLab
Nitrogen	Phthalate free with flow to 30 lpm @ sea level
Maximum Outlet Pressure	116 psig (8 barg)
Hydrocarbon Content	< 2ppm (excluding methane)
Atmospheric Dewpoint	-58°F (-50°C)
Outlet Port	Female 1/4" NPT
Min/Max Ambient Temperature	50°F/95°F (10°C/35°C)
Electrical Requirements	120Vac/60Hz/20Amp / NEMA 5 - 20 Straight Blade
Dimensions	27.6" h x 12.2" w x 35.4" d (70.1cm x 31cm x 90cm)
Shipping Weight	204 lbs. (92.5 kg)

Here's what your colleagues say:

"We've used the Parker Balston Nitroflow® (combined compressor and nitrogen generator) on our LCMS for 3 years. In just over two years, it more than paid for itself in nitrogen savings, but the real advantages of the nitrogen generator are the continuous supply of high quality nitrogen and the tremendous amount of time saved from not having to check, order and switch high pressure liquid nitrogen tanks."

*Karl J. Dria, PhD.
Assistant Research Scientist
Department of Chemistry and
Chemical Biology
Indiana University-Purdue University
Indianapolis*

SOURCE LC/MS TriGas Generator Series

Model LCMS-5000NA

The Parker Balston SOURCE LCMS-5000NA TriGas Generator is a completely engineered system designed to deliver pure nitrogen for curtain gas, pure zero grade air as gas-1/gas-2 source gases and dry -40°F dew point air as source exhaust. The system is designed to produce gases which meet and exceed the requirements of any LC/MS requiring three independent gases.

The system consists of six functional technologies: Coalescing pre-filtration with timed solenoid drains, self regenerating compressed air dehydration membranes, a proprietary heated catalysis module, elegant self-regenerating nitrogen retentate membranes, high capacity - high sensitivity carbon absorption modules and carefully matched final filtration membrane media. These technologies are integrated to a reliable scroll compressor. 3 year compressor pump warranty included. Carries CE Marking and is compliant to WEEE standards.

The Parker Balston SOURCE LCMS-5000NA TriGas Generator will provide enough gas for a single LC/MS instrument on a continual basis and will completely eliminate dependence, expense and hassle with using high pressure nitrogen and zero air cylinders.

The generator can be connected easily, be located in the lab, and features independent stainless steel output gas ports carefully matched to the instrument. Gas distribution, pressure and flow control are integral to each TriGas generator and therefore requirements for secondary gas pressure or gas management systems are eliminated.

There is no longer any need to use valuable laboratory floor space to store excess cylinder reserves to protect from running out of gas, late or missed cylinder deliveries, transportation interruptions or periods of tight supply. With a Parker Balston SOURCE LCMS-5000NA TriGas Generator, you control all your LC/MS gas supplies.



Model LCMS-5000NA



Features and Benefits

- Generates pure nitrogen, zero air and source exhaust air from compressed air, with 24/7 operation
- Eliminates costly and inconvenient nitrogen gas and zero air gas cylinders; minimal annual maintenance
- Prevents running out of gases during LC/MS instrument operation
- Floor standing unit preserves valuable laboratory space and maximizes LC/MS instrument uptime
- Reliable scroll compressor, quiet 49 dB(A) operation at a safe, low pressure
- Gas purity to 99.999% and no phthalates
- Turnkey system that eliminates stainless steel regulators and gas distribution panels
- Produced and manufactured by an ISO 9001 registered organization

Here's what your colleagues say:

"Using the Parker Balston Tri-Gas system with our new ABI 5500 assures maximum uptime of the instrument while offering us the lowest cost to supply gas...The more samples we run lowers our costs and shortens our overall return on the instrument."

*Ed Dabrea
Laboratory Director
Jupiter Environmental Laboratories*

SOURCE LC/MS TriGas Generator Series

Model LCMS-5000NA

Principal Specifications

Model	LCMS-5000NA
Curtain gas (nitrogen)	to 10 lpm and 80 psi
Source gas (uhp zero grade air)	to 23 lpm and 110 psi
Exhaust gas (dry air)	to 8 lpm and 60 psi
Compressor included	Yes - Scroll
Atmospheric dewpoint	-40°F
Hydrocarbons	<0.1 ppm measured as methane
Particles > 0.01 micron	None
Phthalates	None
Suspended liquids	None
Outlets	1/4" tube - stainless steel - 3 each
Dimensions	34"D x 41"W x 43"H
Pressure gauges	3 each
Electrical requirements (1)	120vac, 60Hz, 15 amp and 220vac, 60Hz, 30 amp
Noise level	< 49 dB(A)
Weight	611 lbs. (277 kgs)

NOTES

1 Refer to voltage appendix for electrical and plug configurations for outside North America.

Ordering Information

for assistance, call 800-343-4048, 8 to 5 Eastern Time

Description	Model
Source LC/MS Trigas Generator	LCMS-5000NA
Backup Gas Cylinder Panel	LCMS-EZLINK
Installation Kit	IKLCMS-5000
Preventive Maintenance Plan	LCMS-5000NA-PM
Voltage Reducing Transformer	A03-0286
3/8" Clear PFA Tubing for Remote Compressor Use	11425-1 (specify length)
1/4" Clear PFA Tubing for Connections to LC/MS	11426-1 (specify length)

SOURCE LC/MS TriGas Generator Series Model LCMS-5001NTNA

The Parker Balston SOURCE LCMS TriGas Generator is a completely engineered system designed to transform ordinary compressed air into pure nitrogen for curtain gas, pure zero grade air as gas-1/gas-2 source gases and dry -40°F dew point air as source exhaust. The system is designed to produce gases which meet and exceed the requirements of any LC/MS requiring three independent gases. The system consists of six functional technologies: Coalescing pre-filtration with timed solenoid drains, self regenerating compressed air dehydration membranes, a proprietary heated catalysis module, elegant self-regenerating nitrogen retentate membranes, high capacity - high sensitivity carbon absorption modules and carefully matched final filtration membrane media.

The Parker Balston SOURCE LCMS TriGas Generator will provide enough gas for a single LC/MS instrument on a continual

basis and will completely eliminate dependence, expense and hassle with using high pressure nitrogen and zero air cylinders.

The generator can connect easily to an existing compressed air supply line and features independent stainless steel output gas ports carefully matched to the instrument. Gas distribution, pressure and flow control are integral to each TriGas generator and therefore requirements for secondary gas pressure or gas management systems are eliminated.

There is no longer any need to use valuable laboratory floor space to store excess cylinder reserves to protect from running out of gas, late or missed cylinder deliveries, transportation interruptions or periods of tight supply. With a Parker Balston SOURCE LCMS TriGas Generator, you control all your LC/MS gas supplies.



Model LCMS-5001NTNA



Features and Benefits

- Generates pure nitrogen, zero air and source exhaust air from compressed air; continuously, 24/7
- Eliminates costly and inconvenient nitrogen gas and zero air gas cylinders
- Prevents running out of gases during LC/MS instrument operation
- Floor standing unit preserves valuable laboratory space and maximizes LC/MS instrument uptime
- Reliable, silent operation at a safe, low pressure; minimal annual maintenance
- Gas purity to 99.999% and no Phthalates
- Turnkey system that eliminates stainless steel regulators and gas distribution panels
- Produced and manufactured by an ISO 9001 registered organization
- Listed to U.S. & Canadian safety standards
- Carries CE Marking/compliant to WEEE standard



SOURCE LC/MS TriGas Generator Series

Model LCMS-5001NTNA

Principal Specifications

Model	LCMS-5001NTNA
Curtain gas (nitrogen)	to 10 lpm and 80 psi
Source gas (uhp zero grade air)	to 23 lpm and 110 psi
Exhaust gas (dry air)	to 8 lpm and 60 psi
Air pressure required	85-145 psi (> 100 psi suggested)
Pressure dewpoint	-40°F
Hydrocarbons	<0.1 ppm measured as methane
Particles > 0.01 micron	None
Phthalates	None
Suspended Liquids	None
Inlet	3/8" tube (presto)
Outlets	1/4" tube - stainless steel - 3 each
Dimensions	21"D x 23"W x 41"H
Pressure gauges	3 each
Electrical requirements (1)	120vac, 60Hz, 3 amp
Noise	Silent operation
Weight	157 lbs. (71 kg)

NOTES

1 Refer to voltage appendix for electrical and plug configurations for outside North America.

Ordering Information

for assistance, call 800-343-4048, 8 to 5 Eastern Time

Description	Model
Source LC/MS Trigas Generator	LCMS-5001NTNA
Backup Gas Cylinder Panel	LCMS-EZLINK
Installation Kit	IKLCMS-5000
Preventive Maintenance Plan	LCMS-5001NTNA-PM
Extended Support with 24 Month Warranty	LCMS-5001NT-DN2
1/4" Clear PFA Tubing for Connections to LC/MS	11426-1 (specify length)

SOURCE LC/MS TriGas Generator Series

Model LCMS-5001TNA

The Parker Balston SOURCE LCMS TriGas Generator is a completely engineered system designed to transform ordinary compressed air into pure nitrogen for curtain gas, pure zero grade air as gas-1/gas-2 source gases and dry -40°F dew point air as source exhaust. The system is designed to produce gases which meet and exceed the requirements of any LC/MS requiring three independent gases. The system consists of six functional technologies: Coalescing pre-filtration with timed solenoid drains, self regenerating compressed air dehydration membranes, a proprietary heated catalysis module, elegant self-regenerating nitrogen retentate membranes, high capacity - high sensitivity carbon absorption modules and carefully matched final filtration membrane media.

The Parker Balston SOURCE LCMS TriGas Generator will provide enough gas for a single

LC/MS instrument on a continual basis and will completely eliminate dependence, expense and hassle with using high pressure nitrogen and zero air cylinders.

The generator can connect easily to an existing compressed air supply line and features independent stainless steel output gas ports carefully matched to the instrument. Gas distribution, pressure and flow control are integral to each TriGas generator and therefore requirements for secondary gas pressure or gas management systems are eliminated.

There is no longer any need to use valuable laboratory floor space to store excess cylinder reserves to protect from running out of gas, late or missed cylinder deliveries, transportation interruptions or periods of tight supply. With a Parker Balston SOURCE LCMS TriGas Generator, you control all your LC/MS gas supplies.



Model LCMS-5001NA



Features and Benefits

- Generates pure nitrogen, zero air and source exhaust air from compressed air continuously, 24/7
- Eliminates costly and inconvenient nitrogen gas and zero air gas cylinders
- Prevents running out of gases during LC/MS instrument operation
- Preserves valuable laboratory space and maximizes LC/MS instrument uptime
- Easy installation; reliable, silent operation at a safe, low pressure
- Gas purity to 99.999% and no Phthalates
- Turnkey system that eliminates stainless steel regulators and gas distribution panels
- Produced and manufactured by an ISO 9001 registered organization
- Floor standing, includes internal economizer air receiver system
- Listed to U.S. & Canadian safety standards
- Carries CE Marking/Compliant to WEEE standard

SOURCE LC/MS TriGas Generator Series

Model LCMS-5001TNA

Principal Specifications

Model	LCMS-5001TNA
Curtain gas (nitrogen)	to 10 lpm and 80 psi
Source gas (uhp zero grade air)	to 23 lpm and 110 psi
Exhaust gas (dry air)	to 8 lpm and 60 psi
Air pressure required	85-145 psi (>100 psi suggested)
Pressure dewpoint	-40°F
Hydrocarbons	<0.1 ppm measured as methane
Particles > 0.01 micron	None
Phthalates	None
Suspended Liquids	None
Inlet	3/8" tubing (presto)
Outlets	1/4" tube - stainless steel - 3 each
Dimensions	25"D x 20"W x 43"H
Pressure gauges	3 each
Electrical requirements (1)	120vac, 60Hz, 3 amp
Noise	Silent operation
Weight	271 lbs (123 kgs)

NOTES

1 Refer to voltage appendix for electrical and plug configurations for outside North America.

Ordering Information

for assistance, call 800-343-4048, 8 to 5 Eastern Time

Description	Model
Source Trigas Generator	LCMS-5001TNA
Backup Gas Cylinder Panel	LCMS-EZLINK
Installation Kit	IKLCMS-5000
Preventive Maintenance Plan	LCMS-5001TNA-PM
Extended Support with 24 Month Warranty	LCMS-5001T-DN2
1/4" Clear PFA Tubing for Connections to LC/MS	11426-1 (specify length)

SOURCE LC/MS Super Flow TriGas Generator Series Model LCMS-SF5000NA

The Parker Balston SOURCE LCMS-SF5000NA Super Flow TriGas Generator is a completely engineered system designed to deliver pure nitrogen for curtain gas, pure zero grade air as gas-1/ gas-2 source gases and dry -40°F dew point air as source exhaust. The system is designed to produce gases which meet and exceed the requirements of any LC/MS requiring three independent gases. The system consists of six functional technologies: Coalescing pre-filtration with timed solenoid drains, self regenerating compressed air dehydration membranes, a proprietary heated catalysis module, elegant self-regenerating nitrogen retentate membranes, high capacity - high sensitivity carbon absorption modules and carefully matched final filtration membrane media. These technologies are integrated to a reliable scroll compressor.

The Parker Balston SOURCE LCMS-SF5000NA Super Flow TriGas Generator will provide

enough gas for a single LC/MS instrument on a continual basis and will completely eliminate dependence, expense and hassle with using high pressure nitrogen and zero air cylinders. The generator can be connected easily, be located in the lab, and features independent stainless steel output gas ports carefully matched to the instrument. Gas distribution, pressure and flow control are integral to each TriGas generator and therefore requirements for secondary gas pressure or gas management systems are eliminated.

There is no longer any need to use valuable laboratory floor space to store excess cylinder reserves to protect from running out of gas, late or missed cylinder deliveries, transportation interruptions or periods of tight supply. With a Parker Balston SOURCE LCMS-SF5000NA Super Flow TriGas Generator, you control all your LC/MS gas supplies.



Model LCMS-SF5000NA



Features and Benefits

- Generates pure nitrogen, zero air and source exhaust air from compressed air with continuous 24/7 operation
- 3 year compressor pump warranty
- Prevents running out of gases during LC/MS instrument operation
- Floor standing unit preserves valuable laboratory space and maximizes LC/MS instrument uptime
- Reliable scroll compressor, quiet 49 dB(A) operation at a safe, low pressure
- Gas purity to 99.999% and no phthalates
- Turnkey system that eliminates stainless steel regulators and gas distribution panels; easy installation and minimal maintenance
- Produced and manufactured by an ISO 9001 registered organization
- Carries CE Marking/compliant to WEEE standard

SOURCE LC/MS Super Flow TriGas Generator Series Model LCMS-SF5000NA

Principal Specifications

Model	LCMS-SF5000NA
Curtain gas (nitrogen)	to 20 lpm and 80 psi
Source gas (uhp zero grade air)	to 46 lpm and 110 psi
Exhaust gas (dry air)	to 16 lpm and 60 psi
Compressor included	Yes - Scroll
Atmospheric dewpoint	-40°F
Hydrocarbons	<0.1 ppm measured as methane
Particles > 0.01 micron	None
Phthalates	None
Suspended liquids	None
Outlets	1/4" tube - stainless steel - 6 each
Dimensions	34"D x 61"W x 43"H
Pressure gauges	6 each
Electrical requirements (1)	120vac, 60Hz, 15 amp and 220vac, 60Hz, 30 amp
Noise level	< 49 dB(A)
Weight	788 lbs. (357 kg)

NOTES

1 Refer to voltage appendix for electrical and plug configurations for outside North America.

Ordering Information

for assistance, call 800-343-4048, 8 to 5 Eastern Time

Description	Model
Source LC/MS Trigas Generator	LCMS-SF5000NA
Backup Gas Cylinder Panel	LCMS-SFEZLINK
Installation Kit	IKLCMS-5000
Preventive Maintenance Plan	LCMS-SF5000NA-PM
Voltage Reducing Transformer	A03-0286
3/8" Clear PFA for Remote Compressor Use	11425-1 (specify length)
1/4" Clear PFA Tubing for Connections to LC/MS	11426-1 (specify length)

Low and Mid Flow Nitrogen Generators

Parker Balston® Low Flow Nitrogen Generators include models N2-04, N2-14, N2-14A that produce up to 61 SLPM of compressed nitrogen, on-site. The Parker Balston® Mid-Flow Nitrogen Generators include models N2-22, N2-22ANA, N2-35, and N2-35ANA that produce 132 SLPM of compressed nitrogen, on-site. The purity level of the nitrogen stream is defined by the user, for the application, and may range from 95% to 99.5%.

Low Flow Model N2-14ANA and Mid Flow Models N2-22ANA and N2-35ANA Nitrogen Generators include an oxygen analyzer which monitors the oxygen concentration of the nitrogen stream. An audible alarm signals high or low oxygen concentrations. Parker Balston Nitrogen Generators are complete systems engineered to transform standard compressed air into nitrogen at safe, regulated pressures, on

demand, without the need for operator attention. The systems eliminate the need for costly, dangerous dewars and cylinders in the laboratory.

Nitrogen is produced by utilizing a combination of filtration and membrane separation technologies. A high efficiency prefiltration system pretreats the compressed air to remove all contaminants down to 0.01 micron. Hollow fiber membranes subsequently separate the clean air into a concentrated nitrogen output stream and an oxygen enriched permeate stream, which is vented from the system. The combination of these technologies produces a continuous on demand supply of pure nitrogen.

Typical applications include: LC/MS, nebulizer gas, chemical and solvent evaporation, instrument purge and supply, evaporative light scattering detector use (ELSD), and sparging.



Model N2-22 Mid Flow Membrane Nitrogen Generator




Features and Benefits

- Recommended and used by all major LC/MS manufacturers
- Eliminates the need for costly, dangerous, inconvenient nitrogen cylinders in the laboratory
- Models N2-04, N2-14, N2-22, N2-35 require no electricity
- Compact design frees up valuable laboratory floor space
- Phthalate-free, no organic vapors
- Unlike PSA technology, membrane will not suppress corona needle discharge.



This Technology Features Advanced HiFluxx Fiber



Nitrogen Purity / Flow Chart

Flow measured in SLPM at indicated Operating Pressure, psig. Flows for Model N2-04 printed in black, flows for Models N2-14 and N2-14A in red.

	145	125	110	100	90	80	70	60
99.5	– 11	– 10	– 9	– 8	– 7	– 6	– 5	– 4
99	6 18	5 16	5 15	4 13	4 11	3 10	3 8	2 7
98	11 29	10 25	9 25	8 20	7 18	6 16	5 13	4 11
97	15 40	13 34	13 33	10 27	9 25	8 21	7 18	6 15
96	20 50	17 43	16 42	13 34	12 31	10 26	9 22	7 19
95	24 60	21 52	20 51	17 42	15 37	13 32	11 28	9 24

Nitrogen Purity / Flow Chart

Flow measured in SLPM at indicated Operating Pressure, psig. Flows for Model N2-22, N2-22A printed in black, flows for Models N2-35, N2-35A in red.

	145	125	110	100	90	80	70	60
99.5	19 29	16 25	14 22	13 20	12 18	10 16	9 13	17 11
99	29 44	25 37	22 33	20 30	18 27	15 23	13 20	11 17
98	44 66	38 57	34 51	30 46	27 41	24 36	20 30	17 26
97	59 83	50 74	45 65	40 57	36 52	31 46	26 40	23 35
96	73 109	63 94	56 84	50 75	45 67	39 59	32 50	27 42
95	88 131	77 114	69 102	61 90	55 81	48 71	41 60	35 52

Principal Specifications

Models	N2-04, N2-14, N2-14ANA, N2-22, N2-22ANA, N2-35 and N2-35ANA
Nitrogen Purity	95.0% - 99.5%
Atmospheric Dewpoint	-58°F (-50°C)
Suspended Liquids	None
Particles > 0.01µm	None
Commercially Sterile	Yes
Phthalate-free	Yes
Hydrocarbon-free	Yes
Min./Max. Operating Pressure	60/145 psig
Max. Press. Drop @ 99% N ₂ Purity, 125 psig	10 psig
Recommended Ambient Operating Temperature	68°F (20°C)
Max. Inlet Air Temperature	110°F (43°C)
Inlet/Outlet Ports	1/4" NPT
Electrical Requirements	N2-04, N2-14, N2-22, N2-35 N2-14ANA, N2-22ANA, N2-35ANA None 120 VAC/60 Hz/25 Watts (1)
Shipping Weight	N2-04 N2-14 N2-14ANA, N2-22, N2-22ANA N2-35, N2-35ANA 42.5 lbs (19 kg) 75 lbs (34 kg) 80 lbs (36 kg) 90 lbs (41 kg)
Oxygen Analyzer	Included with Model N2-14ANA, N2-22ANA, N2-35ANA
Dimensions, N2-04	16.1"h x 10.7"w x 13.4"d (40.9cm x 27.2cm x 34cm)
Dimensions, N2-14, N2-14ANA, N2-22, N2-22ANA, N2-35, N2-35ANA	51.5"h x 18"w x 16.2"d (130.8cm x 45.7cm x 41.1cm)

NOTES

1 Refer to voltage appendix for electrical and plug configurations for outside North America.

Ordering Information

for assistance, call 800-343-4048, 8 to 5 Eastern Time

Description	Galvanic Cell	Annual Maintenance Kit	Installation Kit	Preventive Maintenance Plan	Extended Support with 24 Month Warranty
N2-04	N/A	MK7840	IK7572	N2-04 -PM	N2-04-DN2
N2-14	N/A	MK7572C	IK7572	N2-14-PM	N2-14-DN2
N2-14ANA	72695A	MK7572C	IK7572	N2-14A-PM	N2-14A-DN2
N2-22, N2-35	N/A	MK7572C	IK7572	N2-22-PM, N2-35-PM	N2-22-DN2, N2-35-DN2
N2-22ANA, N2-35ANA	72695A	MK7572C	IK7572	N2-22A-PM, N2-35A-PM	N2-22A-DN2, N2-35A-DN2

High Flow Nitrogen Generators

Parker Balston® High Flow Nitrogen Generators include models N2-45, N2-80, N2-135 that produce up to 467 SLPM of compressed nitrogen, on-site. The purity level of the nitrogen stream is defined by the user, for the application, and may range from 95% to 99.5%.

High Flow Model N2-45ANA, N2-80ANA, and N2 135ANA Nitrogen Generators include an oxygen analyzer which monitors the oxygen concentration of the nitrogen stream. An audible alarm signals high or low oxygen concentrations. Parker Balston Nitrogen Generators are complete systems engineered to transform standard compressed air into nitrogen at safe, regulated pressures, on demand, without the need for operator attention. The systems eliminate the need for costly, dangerous dewars and cylinders in the laboratory.

Nitrogen is produced by utilizing a combination of filtration and membrane separation technologies. A high efficiency prefiltration system pretreats the compressed air to remove all contaminants down to 0.01 micron. Hollow fiber membranes subsequently separate the clean air into a concentrated nitrogen output stream and an oxygen enriched permeate stream, which is vented from the system. The combination of these technologies produces a continuous on demand supply of pure nitrogen.

Typical applications include: LC/MS, nebulizer gas, chemical and solvent evaporation, instrument purge and supply, evaporative light scattering detector use (ELSD), and sparging.



Model N2-135 High Flow Membrane Nitrogen Generator



Features and Benefits

- Recommended and used by all major LC/MS manufacturers
- Eliminates the need for costly, dangerous, inconvenient nitrogen cylinders in the laboratory
- Models N2-45, N2-80, and N2-135 require no electricity
- Compact design frees up valuable laboratory floor space
- Phthalate-free, no organic vapors
- Unlike PSA technology, membrane will not suppress corona needle discharge.



High Flow Nitrogen Generators

Nitrogen Purity / Flow Chart

Flow LPM (liters per minute), at 68°F (25°C) inlet air temperature and operating pressure, PSIG.

Flows printed in black are for Models N2-45 and N2-45A

Flows printed in red are for Models N2-80 and N2-80A

Flows printed in green are for Models N2-135 and N2-135A

	145		125		110		100		90		80							
99.5	67	100	133	55	83	110	47	71	94	39	59	78	33	50	66	27	41	54
99	92	138	183	74	112	149	63	95	127	53	79	106	44	66	89	35	53	71
98	129	194	258	106	159	212	89	134	179	73	110	147	62	93	124	50	75	101
97	163	244	325	132	198	264	113	169	226	94	141	187	79	119	159	65	97	130
96	200	300	400	160	240	320	137	205	274	114	171	228	97	145	194	80	119	159
95	233	350	467	187	281	374	160	241	321	134	201	268	111	167	222	90	135	180

Principal Specifications

Model	N2-45, N2-80, N2-135, N2-45ANA, N2-80ANA, and N2-135ANA	
Nitrogen Purity	95.0% - 99.5%	
Atmospheric Dewpoint	-58°F (-50°C)	
Suspended Liquids	None	
Particles > 0.01µm	None	
Commercially Sterile	Yes	
Phthalate-free	Yes	
Hydrocarbon-free	Yes	
Min./Max. Operating Pressure	60/145 psig	
Max. Press. Drop @ 99% N ₂ Purity, 125 psig	10 psig	
Recommended Ambient Operating Temperature	72°F (22°C)	
Max. Inlet Air Temperature	110°F (43°C)	
Inlet/Outlet Ports	1/2" NPT	
Electrical Requirements	N2-45, N2-80, N2-135 N2-45ANA, N2-80ANA, N2-135ANA	None 120 VAC/60 Hz/25 Watts (1)
Shipping Weight	N2-45, N2-80, N2-135 N2-45ANA, N2-80ANA, N2-135ANA	250 lbs (114 kg) 250 lbs (114 kg)
Oxygen Analyzer	Included with Model N2-45ANA, N2-80ANA, N2-135ANA	
Dimensions	67"h x 24"w x 20"d (140cm x 61cm x 50cm)	

NOTES

1 Refer to voltage appendix for electrical and plug configurations for outside North America.

Ordering Information

for assistance, call 800-343-4048, 8 to 5 Eastern Time

Description	Galvanic Cell	Carbon Tower	Maintenance Kit	Installation Kit	Preventive Maintenance Plan	Extended Support with 24 Month Warranty
N2-45	N/A	75344	75478	IK75880	N2-45-PM	N2-45-DN2
N2-45ANA	72695A	75344	75478	IK75880	N2-45A-PM	N2-45A-DN2
N2-80	N/A	75344	75478	IK75880	N2-80-PM	N2-80-DN2
N2-80ANA	72695A	75344	75478	IK75880	N2-80A-PM	N2-80A-DN2
N2-135	N/A	75344	75478	IK75880	N2-135-PM	N2-135-DN2
N2-135ANA	72695A	75344	75478	IK75880	N2-135A-PM	N2-135A-DN2

High Flow Nitrogen Generators

Parker Balston® Monobed Nitrogen Generators produce up to 99.95% pure, compressed nitrogen at dewpoints to -70°F (-21°C) from nearly any compressed air supply. The generators are designed to continually transform standard compressed air into nitrogen at safe, regulated pressures without operator attention.

Parker Balston PSA Nitrogen Generators utilize a combination of filtration and pressure swing adsorption technologies. High efficiency prefiltration pretreats the compressed air to remove all contaminants down to 0.1 micron. Air entering the generator consists of 21% oxygen and 78% nitrogen. The gas separation process preferentially adsorbs oxygen over nitrogen using carbon molecular sieve (CMS). At high pressures the CMS has a greater affinity for oxygen, carbon dioxide, and water vapor than it does at low pressures. By raising and lowering the pressure within the CMS bed, all contami-

nants are captured and released, leaving the CMS unchanged. This process allows the nitrogen to pass through as a product gas at pressure. The depressurization phase of the CMS releases the absorbed oxygen and other contaminant gases to the atmosphere.

The Parker Balston PSA Nitrogen Generators completely eliminate the inconvenience and the high costs of nitrogen Dewars and cylinders. There is no need to depend on outside vendors for your nitrogen gas supplies. The hassles of changing dangerous, high pressure cylinders, and interruption of gas supplies are completely eliminated. The Balston PSA Nitrogen Generators offer long term cost stability eliminating uncontrollable vendor price increases, contract negotiations, long term commitments, and tank rentals. Once the generator is installed, a continuous nitrogen supply of consistent purity is available within minutes from start-up.




Parker Balston
Dual Bed Nitrogen Generators



Features and Benefits

- Lower cost...eliminates the need for costly gas cylinders
- Complete package with prefilters, final filters, and receiving tank
- Compact - frees up valuable floor space
- Eliminates unexpected shutdowns due to a “bad” or empty cylinder
- Hassle-free, easy to install, easy to operate
- Safe and reliable



Pressure swing adsorption gas separation process adsorbs oxygen over nitrogen using carbon molecular sieve (CMS).

High Flow Nitrogen Generators

Parker Balston® Monobed Nitrogen Generators produce up to 99.95% pure, compressed nitrogen at dewpoints to -70°F (-21°C) from nearly any compressed air supply. The generators are designed to continually transform standard compressed air into nitrogen at safe, regulated pressures without operator attention.

Parker Balston PSA Nitrogen Generators utilize a combination of filtration and pressure swing adsorption technologies. High efficiency prefiltration pretreats the compressed air to remove all contaminants down to 0.1 micron. Air entering the generator consists of 21% oxygen and 78% nitrogen. The gas separation process preferentially adsorbs oxygen over nitrogen using carbon molecular sieve (CMS). At high pressures the CMS has a greater affinity for oxygen, carbon dioxide, and water vapor than it does at low pressures.

By raising and lowering the pressure within the CMS bed, all contaminants are captured and released, leaving the CMS unchanged. This process allows the nitrogen to pass through as a product gas at pressure. The depressurization phase of the CMS releases the absorbed oxygen and other contaminant gases to the atmosphere.

The Parker Balston PSA Nitrogen Generators completely eliminate the inconvenience and the high costs of nitrogen Dewars and cylinders. There is no need to depend on outside vendors for your nitrogen gas supplies. The hassles of changing dangerous, high pressure cylinders, and interruption of gas supplies are completely eliminated. The Balston PSA Nitrogen Generators offer long term cost stability eliminating uncontrollable vendor price increases, contract negotiations, long term commitments, and tank rentals. Once the generator is installed, a continuous nitrogen supply of consistent purity is available within minutes from start-up.

Installation consists of simply connecting a standard compressed air line to the inlet and connecting the outlet to a nitrogen line. Plug the electrical cord into a wall outlet, and the unit is ready for trouble-free operation. This system is designed to operate 24 hours per day, 7 days per week.

Once the system is operating, it requires little monitoring. The only maintenance involves changing the coalescing prefilter cartridges and final sterile air filter periodically. The PSA towers do not require any maintenance.

An oxygen monitor to measure the oxygen concentration of the nitrogen stream is available as an option. An audible alarm signals high or low oxygen concentrations (determined by the application). The oxygen analyzer is supplied with alarm relay outputs which may be used to signal a remote alarm, open a backup supply or the process stream, or close the process flow for protection of downstream equipment or processes.

Principal Specifications

Model	AGS200, AGS400	
Nominal Conditions		
Feed Pressure	140 psig	
Temperature	80°F	
Ambient Pressure	1 Atm.	
Compressed Air Specifications		
Maximum Pressure	140 psig	
Temperature Range	60°F - 105°F	
Dewpoint	40°F pressure dewpoint or better	
Residual Oil Content	Trace	
Particles	<.01 micron	
Ambient Conditions		
Temperature	45°F-90°F	
Ambient Pressure	Atmospheric	
Air Quality	Clean air without contaminants	
Dimensions	28.5"L x 32.25"D x 76.25"H	
Weight	AGS200 AGS400	520 lbs (236 kg) 738 lbs (335 kg)
Inlet	1/2" NPT	
Outlet	1/2" NPT	

Nitrogen Purity Flow Chart

	Flow Rate (SCFH)	Flow Rate (SCFH)
Model	99.9%, 140 psig	99.99%, 140 psig
AGS200	235	47
AGS400	470	94

High Flow Nitrogen Generators

Parker Balston® High Flow Nitrosource Nitrogen Generators produce up to 99.5% pure, commercially sterile nitrogen at dewpoints to -58°F (-50°C) from a compressed air supply. All Membrane Nitrogen Generators include a 0.01 micron membrane filter which ensures the nitrogen is completely free of suspended impurities.

Parker Balston High Flow Nitrosource Nitrogen Generators are one of the most efficient membrane systems available with higher recovery rates and lower operating costs than many other membrane systems.

The generators utilize proprietary membrane separation technology. The membrane divides the air into two separate streams: one is 95%-99.5% pure nitrogen, and the other is oxygen rich with carbon dioxide and other trace gases.

The generator separates air into its component gases by passing inex-

ensive, conventional compressed air through bundles of individual hollow fiber, semi-permeable membranes. Each fiber has a perfectly circular cross section and a uniform bore through its center. Because the fibers are so small, a great many can be packed into a limited space, providing an extremely large membrane surface area that can produce a relatively high volume product stream.

Compressed air is introduced to the center of the fibers at one end of the module and contacts the membrane as it flows through the fiber bores.

While oxygen, water vapor and other trace gases permeate the membrane fiber and are discharged through a permeate port, the nitrogen is contained within the hollow fiber membrane, and flows through the outlet port of the module.

Water vapor also permeates through the membrane; therefore, the nitrogen product gas is very dry.



Parker Balston N2-300
Nitrosource Nitrogen Generators



Features and Benefits

- Lower cost...eliminates the need for costly gas cylinders
- Complete package with prefilters, carbon filter, and membrane filter
- Compact - frees up valuable floor space
- Eliminates unexpected shutdowns due to a "bad" or empty cylinder
- Hassle-free, easy to install, easy to operate
- Safe and reliable
- Expandable modular design

Applications

High thru-put LC/MS contract labs
Sample concentrators
Nitrogen supply to analytical lab

Custom Systems Available

Flow rates to 2,265 lpm
Delivery pressures to customer's specifications
Skid mounted systems with compressor, receiving tank and controls are available

High Flow Nitrogen Generators

The Parker Balston Nitrosource Nitrogen Generators completely eliminate and inconvenient and the high costs of nitrogen Dewars and cylinders. There is no need to depend on outside vendors for nitrogen gas supplies. The hassles of changing dangerous, high pressure cylinders and interruption of gas supplies are completely eliminated. The Balston Systems offer long term cost stability by eliminating uncontrollable vendor price increases, contract negotiation, long term commitments and tank rentals. Once the generator is installed, a continuous nitrogen supply of consistent purity is available within minutes from start-up.

The Parker Balston Nitrosource Nitrogen Generators are com-

plete systems ready to operate as delivered with carefully matched components engineered for easy installation, operation and long term reliability.

The generators are free-standing and housed in an attractive cabinet. Standard features include: high efficiency coalescing prefilters with automatic drains, an activated carbon filter, and a 0.01 micron membrane final filter. Installation consists of simply connecting a standard compressed air line to the inlet and connecting the outlet to a nitrogen line.

There is no complicated operating procedure to learn or labor intensive monitoring involved. Simply select the purity your process requires, set the flow and pressure, and within

minutes high purity, dry nitrogen is available for use!

Once the system is operating, it requires little monitoring. The only maintenance involves changing the coalescing filter cartridges and activated carbon filter periodically. This is a simple ten minute procedure.

All models also include an oxygen monitor which offers LCD readouts and remote alarm or chart recorder capabilities. An audible alarm signals high or low oxygen concentrations (determined by the application). The oxygen monitor is supplied with alarm relay outputs which may be used to signal a remote alarm, open a backup supply or the process stream, or close the process flow.

Principal Specifications - Nitrosource Series

Model	N2-300	N2-460	N2-600
Atmospheric Dewpoint	-58°F (-50°C)	-58°F (-50°C)	-58°F (-50°C)
Commercially Sterile	Yes	Yes	Yes
Particles >0.01 micron	None	None	None
Suspended Liquids	None	None	None
Min/Max Operating Pressure	60 psig/145 psig	60 psig/145 psig	60 psig/145 psig
Max Pressure Drop (at 95% N ₂ , 125 psig)	15 psig	15 psig	15 psig
Operating Temperature	70°F (21°C)	70°F (21°C)	70°F (21°C)
Min/Max Inlet Air Temp.	50°F /104°F (10°F /40°F)	50°F /104°F (10°F /40°F)	50°F /104°F (10°F /40°F)
Recommended Inlet Air Temp.	70°F (21°C)	70°F (21°C)	70°F (21°C)
Electrical Requirements	90-250 VAC 50-60 Hz	90-250 VAC 50-60 Hz	90-250 VAC 50-60 Hz
Dimensions	29"W x 31"D x 76"H (74cm x 51cm x 193cm)	29"W x 42"D x 76"H (74cm x 79cm x 193cm)	29"W x 53"D x 76"H (74cm x 107cm x 193cm)
Shipping Weight	660 lbs. (300 kg)	870 lbs. (395 kg)	1,290 lbs. (586 kg)

Flow Rates (lpm) @ 100 psig, 68°F

Model	99.5%	99%	98%	97%	96%	95%
N2-300	200	311	538	736	935	1133
N2-460	297	467	807	1104	1402	1699
N2-600	396	623	1076	1473	1869	2266

Voltage Appendix

220vac / 50hz configuration for locations where final plug configuration is unknown

Order Part Number



FID-1000-220, FID-2500-220, FID-3500-220, GCGS-7890-220, H2PD-150-220, H2PD-300-220, 75-83-220, HPZA-3500-220, HPZA-7000-220, HPZA18000-220, HPZA30000-220, HPN2-1100-220, HPN2-2000-220, UHPN2-1100-220, 76-97-220, 76-98-220, 74-5041-220, UDA-300-220, LCMS-5000-220, LCMS-5001T-220, LCMS-5001NT-220, N2-14A, N2-22A, N2-35A, N2-45A, N2-80A, N2-135A, MGG-400-220, MGG-2500-220, TOC-625-220, TOC-1250-220

* Units will be supplied only with IEC connector as depicted, power cord to be customer supplied

220vac / 50hz plug configuration for Australia

Order Part Number



FID-1000AU, FID-2500AU, FID-3500AU, GCGS-7890AU, H2PD-150AU, H2PD-300AU, 75-83AU, HPZA-3500AU, HPZA-7000AU, HPZA-18000AU, HPZA-30000AU, HPN2-1100AU, HPN2-2000AU, UHPN2-1100AU, 76-97AU, 76-98AU, 74-5041AU, UDA-300AU, LCMS-5000AU, LCMS-5001TAU, LCMS-5001NTAU, N2-14AAU, N2-22AAU, N2-35AAU, N2-45AAU, N2-80AAU, N2-135AAU, MGG-400AU, MGG-2500AU, TOC-625AU, TOC-1250AU

* Models 75-45AU, 75-52AU and 75-62AU will include universal fit plug and transformer kit.

220vac / 50hz plug configuration for Europe

Order Part Number



FID-1000EU, FID-2500EU, FID-3500EU, GCGS-7890EU, H2PD-150EU, H2PD-300EU, 75-83EU, HPZA-3500EU, HPZA-7000EU, HPZA-18000EU, HPZA-30000EU, HPN2-1100EU, HPN2-2000EU, UHPN2-1100EU, 76-97EU, 76-98EU, 74-5041EU, UDA-300EU, LCMS-5000EU, LCMS-5001TEU, LCMS-5001NTEU, N2-14AEU, N2-22AEU, N2-35AEU, N2-45AEU, N2-80AEU, N2-135AEU, MGG-400EU, MGG-2500EU, TOC-625EU, TOC-1250EU

* Models 75-45EU, 75-52EU and 75-62EU will include universal fit plug and transformer kit.

100vac / 60hz plug configuration for Japan

Order Part Number



FID-1000JA-100, FID-2500JA-100, FID-3500JA-100, GCGS-7890JA-100, H2PD-150JA-100, H2PD-300JA-100, 75-83JA-100, HPZA-3500JA-100, HPZA-7000JA-100, HPZA-18000JA-100, HPZA-30000JA-100, HPN2-1100JA-100, HPN2-2000JA-100, UHPN2-1100JA-100, 76-97JA-100, 76-98JA-100, 74-5041JA-100, UDA-300JA-100, LCMS-5000JA-100, LCMS-5001TJA-100, LCMS-5001NTJA-100, N2-14AJA-100, N2-22AJA-100, N2-35AJA-100, N2-45AJA-100, N2-80AJA-100, N2-135AJA-100, MGG-400JA-100, MGG-2500JA-100, TOC-625JA-100, TOC-1250JA-100

* Models 75-45JA-100, 75-52JA-100 and 75-62JA-100 will include universal fit plug and transformer kit.

220vac / 50hz plug configuration for United Kingdom (some Asia)

Order Part Number



FID-1000UK, FID-2500UK, FID-3500UK, GCGS-7890UK, H2PD-150UK, H2PD-300UK, 75-83UK, HPZA-3500UK, HPZA-7000UK, HPZA-18000UK, HPZA-30000UK, HPN2-1100UK, HPN2-2000UK, UHPN2-1100UK, 76-97UK, 76-98UK, 74-5041UK, UDA-300UK, LCMS-5000UK, LCMS-5001TUK, LCMS-5001NTUK, N2-14AUK, N2-22AUK, N2-35AUK, N2-45AUK, N2-80AUK, N2-135AUK, MGG-400UK, MGG-2500UK, TOC-625UK, TOC-1250UK

* Models 75-45UK, 75-52UK and 75-62UK will include universal fit plug and transformer kit.

Recommended Gas Generators for Analytical Instruments

Instrument	Gas Requirements	Gas Purity Requirements	Flow Rates	Generator Recommendation/Model
Atomic Absorption (AA) with Flame	Air for Oxidant Gas	Clean, Dry	1-7 SCFM	AA Gas Purifier (Model 73-100)
Atomic Thermal Desorber	Zero Air	Clean, Dry, Hydrocarbon-free	Up to 1600 ml/min.	Zero Air or TOC Gas Generator (HPZA-3500 or TOC-1250)
	Hydrogen for FID Fuel	Clean, Dry, High Purity	Up to 40 ml/min. per FID	Hydrogen Generator (H2PEM-100, H2PEM-165) (H2PEM-260, H2PEM-510)
Atmospheric Pressure Ionization (API-MS)	Air for Nebulizer Gas	Clean, Dry, Hydrocarbon-free	< 30 LPM	Zero Air Generator (HPZA-30000)
	Nitrogen for Curtain, Sheath, and Shield gas	99% or higher	< 20 LPM	Nitrogen Generator (N2-14, N2-22, N2-35, NitroFlowLab)
Autosamplers for Various Instruments	Air for Pneumatic Controls	Clean, Dry	< 1 SCFM	Membrane Air Dryer (64-02)
	Nitrogen for Sample Injector	Ultra High Purity	< 550 cc/min	UHP Nitrogen Generator (HPN2-1100) (UHPN2-1100)
CO ₂ Analyzers	Calibration Air	CO ₂ -free	0.5-1.0 SLPM	FT-IR Purge Gas Generator (75-45, 75-52)
Continuous Emissions Monitoring (CEM)	Calibration Air Dilution Air	Dry, CO ₂ , SO ₂ , NO _x , Hydrocarbon-free	10-15 SLPM	CEM Zero Air Generator (75-45-M744)
Emissions Analyzers	Zero Air	Hydrocarbon-free	2-15 SLPM	Zero Air Generator (HPZA-18000)
Fourier Transform Infrared Spectrometer (FT-IR)	Air for Sample Compartment, Optics, and/or Air-Bearing Components	Clean, Dry, CO ₂ -free	0.5-3 SCFM	FT-IR Purge Gas Generator (75-62, 75-52, 75-45) Lab Gas Generator (74-5041NA)
Gas Chromatograph (GC) GC-FID	Zero Air as Flame Support Air	Clean, Hydrocarbon-free	150-600 cc/min.	Zero Air Generator (HPZA-3500)
	Hydrogen as Flame Fuel Gas	Ultra High Purity	30-40 cc/min.	Hydrogen Generator (H2PEM-260)
	Hydrogen as Capillary Carrier Gas	Ultra High Purity	Varies	Hydrogen Generator (H2PD-300)
	Nitrogen as Packed Carrier Gas	Ultra High Purity, Zero Grade	Varies	UHP Nitrogen Generator (UHPN2-1100)
	Nitrogen as Make up Gas	Ultra High Purity, Zero Grade	<100 cc/min	UHP Nitrogen Generator (UHPN2-1100)
GC-FPD	Zero Air as Flame Support Air	Clean, Hydrocarbon-free	<200 cc/min	Zero Air Generator (HPZA-3500)
	Hydrogen as Flame Fuel Gas	Ultra High Purity	50-70 cc/min	Hydrogen Generator (H2PEM-260)
	Hydrogen as Capillary Carrier Gas	Ultra High Purity	Varies	Hydrogen Generator (H2-1200)
	Nitrogen as Packed Carrier Gas	Ultra High Purity	Varies	UHP Nitrogen Generator (UHPN2-1100)
GC-NPD	Zero Air to Rubidium/Thermonic Bead	Dry, Clean, Hydrocarbon-Free	60-200 cc/min.	Zero Air Generator (HPZA-3500)
	Hydrogen as Detector Support Gas	Ultra High Purity	<10 cc/min	Hydrogen Generator (H2PEM-100)
	Hydrogen as Capillary Carrier Gas	Ultra High Purity	Varies	Hydrogen Generator (Palladium) (H2PD-300)
	Nitrogen as Packed Carrier Gas	Ultra High Purity	Varies	UHP Nitrogen Generator (UHPN2-1100)
GC-ECD	Nitrogen as Carrier Gas	Ultra High Purity, Zero Grade	Varies	UHP Nitrogen Generator (UHPN2-1100)
	Nitrogen as Make up Gas	Ultra High Purity, Zero Grade	<100 cc/min	UHP Nitrogen Generator (UHPN2-1100)
GC-ELCD, HALL	Hydrogen as Reaction Gas	Ultra High Purity	70-200 cc/min	Hydrogen Generator (H2PD-300)

Recommended Gas Generators for Analytical Instruments

Instrument	Gas Requirements	Gas Purity Requirements	Flow Rates	Generator Recommendation/Model
GC-TCD	Hydrogen as Carrier & Reference Gas	Ultra High Purity	Varies	Hydrogen Generator (H2PD-300)
LC/MS	Nitrogen as a Curtain Gas	LC/MS Grade	3-30 lpm	Nitrogen Generator (N2-14, N2-14ANA, NitroFlowLab) (NitroFlow60, N2-35, N2-35ANA)
ICP Spectrometer	Nitrogen as Optic/Camera Purge	Ultra High Purity	<1-5 lpm	Nitrogen Generator (76-97NA, 76-98NA)
Nuclear Magnetic Resonance (NMR)	Air for Lifting, Spinning	Clean, Dry	<10 SCFM	Air Dryer (UDA-300NA) Lab Gas Generator (74-5041NA)
Ozone Generator	Supply Air	Clean, Dry	.3-20 SCFM	Air Dryer (64-01, 64-02, 64-10, UDA-300NA)
Protein Analyzer	Dry Air, Nitrogen	Clean, Dry	40 psig	Nitrogen Generator (N2-14, N2-22, NitroFlowLab, N2-35)
Solvent Evaporators (Sample Concentrators)	Nitrogen	Clean, Dry Nitrogen	Up to 5 SCFM	Nitrogen Generator (Nitrovap-1LV, Nitrovap-2LV)
Stack Gas Sampler	Dilution Air	Clean, Dry	<1.0 SCFM	CEM Zero Air Generator (75-45-M744)
Total Oxygen Demand (TOD)	Nitrogen Carrier Gas	Ultra High Purity	300 cc/min	Nitrogen Generator (UHPN2-1100)
Thermal Gravimetric Analyzer (TGA)	Nitrogen as Furnace Purge	Clean, Dry, Inert	<100 cc/min	Nitrogen Generator (UHPN2-1100)
Differential Scanning Calorimeter (DSC)	Air for Air Shield	Clean, Dry	<100 cc/min	Dry Air Generator (64-01, UDA-300)
Total Hydrocarbon Analyzer (THA)	Zero Air for FID	Clean, Hydrocarbon-Free	50-500 cc/min	Zero Air Generator (75-82S, 75-83NA)
	Hydrogen as Flame Fuel Gas	Ultra High Purity	5-50 cc/min	Hydrogen Generator (H2PEM-100)
Total Organic Carbon Analyzer (TOC)	Dry Air or Nitrogen for Carrier Gas or Combustion Gas	Clean, Dry, Hydrocarbon-Free	100-500 SLPM	TOC Gas Generator (TOC-625, TOC-1250)
		CO ₂ -Free Ultra High Purity	50-700 cc/min	UHP Nitrogen Generator (UHPN2-1100)

Parker Balston also offers Gas Generators for these Applications



Products for LC/MS & Evaporation

(Request Bulletin AGS-LCMS)

- High purity nitrogen for LCMS instruments and solvent evaporation
- Tri-gas units available for instruments that require nitrogen, dry air and zero grade air
- Produce a continuous supply of high purity nitrogen from an existing compressed air supply
- Integrated compressor systems eliminate the need for house air
- Systems available to support one or dozens



Products for Chromatography

(Request Bulletin AGS-Chromatography)

- Hydrogen, Zero Air and UHP Nitrogen Generators for Gas Chromatography
- Combination systems available to provide multiple gasses from one unit
- Highest purities available from any supplier



Products for Spectroscopy

(Request Bulletin AGS-Spectroscopy)

- Remove water and CO₂ from compressed air
- Protect expensive optics from damage from water vapor
- Increase Signal to Noise Ratio and maximize instrument sensitivity
- Ultra dry air for NMR injecting, spinning and ejecting samples



Products for TOC Analysis

(Request Bulletin AGS-TOC)

- Generate gasses for all combustion, UV persulfate and wet oxidation techniques
- Ensures consistent, reliable, instrument operation and reduces instrument service and maintenance costs



Products for Ultra Dry Air

(Request Bulletin AGS-JDA)

- Gas generators for dilution and calibration of Emissions Analyzers
- Exceed instrument manufacturer specifications
- Nitrogen and specialty blend gasses available



Analytical Gas Supplies

(Request Bulletin AGS SUPCAT)

- Installation kits, compressors, purifiers, flow-meters, regulators and all the materials needed to equip your lab
- High quality components, designed specifically for use with Parker gas generators, to deliver high purity gas to your instruments

Gas Generator Services



Parker Balston Extended Support Services extend the warranty term of gas generators to 24-months. There are two choices available for level of service: Depot and Express. All parts and labor are included, with “next business morning” delivery available.



Parker Balston “Balston Bucks” Loyalty Programs are offered to every customer who purchases gas generators. Services include special discounts and incentives on gas generator spare parts and consumables as well as special deals on buying your next gas generator. Customers can easily “opt-in” and opt-out” of our special e-mail alerts system which features newsletters, new product announcements and product reminders.



Parker Balston Leasing and Rental Services can provide simple cost effective ways to acquire your next gas generator. Our competitive rates typically provide a monthly payment less than current monthly cylinder gas expenditures. Leasing and rental programs help the customer avoid the need to use capital budget money.



Parker Balston Preventative Maintenance Contracts provide convenient direct in-lab maintenance service for your gas generator. A factory trained technician will service your gas generator, in your lab, with original Parker parts. Preventative maintenance saves time, money and will reduce the total cost of ownership of your gas generator.

Application Notes

Offer of Sale

The items described in this document and other documents and descriptions provided by Parker Hannifin Corporation, its subsidiaries and its authorized distributors ("Seller") are hereby offered for sale at prices to be established by Seller. This offer and its acceptance by any customer ("Buyer") shall be governed by all of the following Terms and Conditions. Buyer's order for any item described in its document, when communicated to Seller verbally, or in writing, shall constitute acceptance of this offer. All goods, services or work described will be referred to as "Products".

- 1. Terms and Conditions.** Seller's willingness to offer Products, or accept an order for Products, to or from Buyer is subject to these Terms and Conditions or any newer version of the terms and conditions found on-line at www.parker.com/saleterms/. Seller objects to any contrary or additional terms or conditions of Buyer's order or any other document issued by Buyer.
- 2. Price Adjustments; Payments.** Prices stated on Seller's quote or other documentation offered by Seller are valid for 30 days, and do not include any sales, use, or other taxes unless specifically stated. Unless otherwise specified by Seller, all prices are F.C.A. Seller's facility (INCOTERMS 2010). Payment is subject to credit approval and is due 30 days from the date of invoice or such other term as required by Seller's Credit Department, after which Buyer shall pay interest on any unpaid invoices at the rate of 1.5% per month or the maximum allowable rate under applicable law.
- 3. Delivery Dates; Title and Risk; Shipment.** All delivery dates are approximate and Seller shall not be responsible for any damages resulting from any delay. Regardless of the manner of shipment, title to any products and risk of loss or damage shall pass to Buyer upon placement of the products with the shipment carrier at Seller's facility. Unless otherwise stated, Seller may exercise its judgment in choosing the carrier and means of delivery. No deferral of shipment at Buyer's request beyond the respective dates indicated will be made except on terms that will indemnify, defend and hold Seller harmless against all loss and additional expense. Buyer shall be responsible for any additional shipping charges incurred by Seller due to Buyer's acts or omissions.
- 4. Warranty.** Seller warrants that the Products sold hereunder shall be free from defects in material or workmanship for a period of twelve months from the date of delivery to Buyer or 2,000 hours of normal use, whichever occurs first. The prices charged for Seller's products are based upon the exclusive limited warranty stated above, and upon the following disclaimer: **DISCLAIMER OF WARRANTY: THIS WARRANTY COMPRISES THE SOLE AND ENTIRE WARRANTY PERTAINING TO PRODUCTS PROVIDED HEREUNDER. SELLER DISCLAIMS ALL OTHER WARRANTIES, EXPRESS AND IMPLIED, INCLUDING DESIGN, MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.**
- 5. Claims; Commencement of Actions.** Buyer shall promptly inspect all Products upon delivery. No claims for shortages will be allowed unless reported to the Seller within 10 days of delivery. No other claims against Seller will be allowed unless asserted in writing within 30 days after delivery. Buyer shall notify Seller of any alleged breach of warranty within 30 days after the date the defect is or should have been discovered by Buyer. Any action based upon breach of this agreement or upon any other claim arising out of this sale (other than an action by Seller for an amount due on any invoice) must be commenced within 12 months from the date of the breach without regard to the date breach is discovered.
- 6. LIMITATION OF LIABILITY.** UPON NOTIFICATION, SELLER WILL, AT ITS OPTION, REPAIR OR REPLACE A DEFECTIVE PRODUCT, OR REFUND THE PURCHASE PRICE. **IN NO EVENT SHALL SELLER BE LIABLE TO BUYER FOR ANY SPECIAL, INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF, OR AS THE RESULT OF, THE SALE, DELIVERY, NON-DELIVERY, SERVICING, USE OR LOSS OF USE OF THE PRODUCTS OR ANY PART THEREOF, OR FOR ANY CHARGES OR EXPENSES OF ANY NATURE INCURRED WITHOUT SELLER'S WRITTEN CONSENT, EVEN IF SELLER HAS BEEN NEGLIGENT, WHETHER IN CONTRACT, TORT OR OTHER LEGAL THEORY. IN NO EVENT SHALL SELLER'S LIABILITY UNDER ANY CLAIM MADE BY BUYER EXCEED THE PURCHASE PRICE OF THE PRODUCTS.**
- 7. User Responsibility.** The user, through its own analysis and testing, is solely responsible for making the final selection of the system and Product and assuring that all performance, endurance, maintenance, safety and warning requirements of the application are met. The user must analyze all aspects of the application and follow applicable industry standards and Product information. If Seller provides Product or system options, the user is responsible for determining that such data and specifications are suitable and sufficient for all applications and reasonably foreseeable uses of the Products or systems.
- 8. Loss to Buyer's Property.** Any designs, tools, patterns, materials, drawings, confidential information or equipment furnished by Buyer or any other items which become Buyer's property, will be considered obsolete and may be destroyed by Seller after two consecutive years have elapsed without Buyer ordering the items manufactured using such property. Seller shall not be responsible for any loss or damage to such property while it is in Seller's possession or control.
- 9. Special Tooling.** A tooling charge may be imposed for any special tooling, including without limitation, dies, fixtures, molds and patterns, acquired to manufacture Products. Such special tooling shall be and remain Seller's property notwithstanding payment of any charges by Buyer. In no event will Buyer acquire any interest in apparatus belonging to Seller which is utilized in the manufacture of the Products, even if such apparatus has been specially converted or adapted for such manufacture and notwithstanding any charges paid by Buyer. Unless otherwise agreed, Seller shall have the right to alter, discard or otherwise dispose of any special tooling or other property in its sole discretion at any time.
- 10. Buyer's Obligation; Rights of Seller.** To secure payment of all sums due or otherwise, Seller shall retain a security interest in the goods delivered and this agreement shall be deemed a Security Agreement under the Uniform Commercial Code. Buyer authorizes Seller as its attorney to execute and file on Buyer's behalf all documents Seller deems necessary to perfect its security interest.
- 11. Improper use and Indemnity.** Buyer shall indemnify, defend, and hold Seller harmless from any claim, liability, damages, lawsuits, and costs (including attorney fees), whether for personal injury, property damage, patent, trademark or copyright infringement or any other claim, brought by or incurred by Buyer, Buyer's employees, or any other person, arising out of: (a) improper selection, improper application or other misuse of Products purchased by Buyer from Seller; (b) any act or omission, negligent or otherwise, of Buyer; (c) Seller's use of patterns, plans, drawings, or specifications furnished by Buyer to manufacture Product; or (d) Buyer's failure to comply with these terms and conditions. Seller shall not indemnify Buyer under any circumstance except as otherwise provided.
- 12. Cancellations and Changes.** Orders shall not be subject to cancellation or change by Buyer for any reason, except with Seller's written consent and upon terms that will indemnify, defend and hold Seller harmless against all direct, incidental and consequential loss or damage. Seller may change product features, specifications, designs and availability with notice to Buyer.
- 13. Limitation on Assignment.** Buyer may not assign its rights or obligations under this agreement without the prior written consent of Seller.
- 14. Force Majeure.** Seller does not assume the risk and shall not be liable for delay or failure to perform any of Seller's obligations by reason of circumstances beyond the reasonable control of Seller (hereinafter "Events of Force Majeure"). Events of Force Majeure shall include without limitation: accidents, strikes or labor disputes, acts of any government or government agency, acts of nature, delays or failures in delivery from carriers or suppliers, shortages of materials, or any other cause beyond Seller's reasonable control.
- 15. Waiver and Severability.** Failure to enforce any provision of this agreement will not waive that provision nor will any such failure prejudice Seller's right to enforce that provision in the future. Invalidation of any provision of this agreement by legislation or other rule of law shall not invalidate any other provision herein. The remaining provisions of this agreement will remain in full force and effect.
- 16. Termination.** Seller may terminate this agreement for any reason and at any time by giving Buyer thirty (30) days written notice of termination. Seller may immediately terminate this agreement, in writing, if Buyer: (a) commits a breach of any provision of this agreement (b) appoints a trustee, receiver or custodian for all or any part of Buyer's property (c) files a petition for relief in bankruptcy on its own behalf, or by a third party (d) makes an assignment for the benefit of creditors, or (e) dissolves or liquidates all or a majority of its assets.
- 17. Governing Law.** This agreement and the sale and delivery of all Products hereunder shall be deemed to have taken place in and shall be governed and construed in accordance with the laws of the State of Ohio, as applicable to contracts executed and wholly performed therein and without regard to conflicts of laws principles. Buyer irrevocably agrees and consents to the exclusive jurisdiction and venue of the courts of Cuyahoga County, Ohio with respect to any dispute, controversy or claim arising out of or relating to this agreement.
- 18. Indemnity for Infringement of Intellectual Property Rights.** Seller shall have no liability for infringement of any patents, trademarks, copyrights, trade dress, trade secrets or similar rights except as provided in this Section. Seller will defend and indemnify Buyer against allegations of infringement of U.S. patents, U.S. trademarks, copyrights, trade dress and trade secrets ("Intellectual Property Rights"). Seller will defend at its expense and will pay the cost of any settlement or damages awarded in an action brought against Buyer based on an allegation that a Product sold pursuant to this Agreement infringes the Intellectual Property Rights of a third party. Seller's obligation to defend and indemnify Buyer is contingent on Buyer notifying Seller within ten (10) days after Buyer becomes aware of such allegations of infringement, and Seller having sole control over the defense of any allegations or actions including all negotiations for settlement or compromise. If a Product is subject to a claim that it infringes the Intellectual Property Rights of a third party, Seller may, at its sole expense and option, procure for Buyer the right to continue using the Product, replace or modify the Product so as to make it noninfringing, or offer to accept return of the Product and return the purchase price less a reasonable allowance for depreciation. Notwithstanding the foregoing, Seller shall have no liability for claims of infringement based on information provided by Buyer, or directed to Products delivered hereunder for which the designs are specified in whole or part by Buyer, or infringements resulting from the modification, combination or use in a system of any Product sold hereunder. The foregoing provisions of this Section shall constitute Seller's sole and exclusive liability and Buyer's sole and exclusive remedy for infringement of Intellectual Property Rights.
- 19. Entire Agreement.** This agreement contains the entire agreement between the Buyer and Seller and constitutes the final, complete and exclusive expression of the terms of sale. All prior or contemporaneous written or oral agreements or negotiations with respect to the subject matter are herein merged.
- 20. Compliance with Law, U. K. Bribery Act and U.S. Foreign Corrupt Practices Act.** Buyer agrees to comply with all applicable laws and regulations, including both those of the United Kingdom and the United States of America, and of the country or countries of the Territory in which Buyer may operate, including without limitation the U. K. Bribery Act, the U.S. Foreign Corrupt Practices Act ("FCPA") and the U.S. Anti-Kickback Act (the "Anti-Kickback Act"), and agrees to indemnify and hold harmless Seller from the consequences of any violation of such provisions by Buyer, its employees or agents. Buyer acknowledges that they are familiar with the provisions of the U. K. Bribery Act, the FCPA and the Anti-Kickback Act, and certifies that Buyer will adhere to the requirements thereof. In particular, Buyer represents and agrees that Buyer shall not make any payment or give anything of value, directly or indirectly to any governmental official, any foreign political party or official thereof, any candidate for foreign political office, or any commercial entity or person, for the purpose of influencing such person to purchase products or otherwise benefit the business of Seller.

Worldwide Filtration Manufacturing Locations

North America

Compressed Air Treatment Filtration & Separation/Balston

Haverhill, MA
978 858 0505
www.parker.com/balston

Finite Airtek Filtration Airtek/domnick hunter/Zander

Lancaster, NY
716 686 6400
www.parker.com/faf

Finite Airtek Filtration/Finite

Oxford, MI
248 628 6400
www.parker.com/finitefilter

Engine Filtration & Water Purification Racor

Modesto, CA
209 521 7860
www.parker.com/racor

Racor

Holly Springs, MS
662 252 2656
www.parker.com/racor

Racor

Beaufort, SC
843 846 3200
www.parker.com/racor

Racor – Village Marine Tec.

Gardena, CA
310 516 9911
desalination.parker.com

Hydraulic Filtration

Hydraulic Filter

Metamora, OH
419 644 4311
www.parker.com/hydraulicfilter

Laval, QC Canada
450 629 9594
www.parkerfarr.com

Process Filtration

domnick hunter Process Filtration

Oxnard, CA
805 604 3400
www.parker.com/processfiltration

Europe

Compressed Air Treatment domnick hunter Filtration & Separation

Gateshead, England
+44 (0) 191 402 9000
www.parker.com/dhfn

Parker Gas Separations

Etten-Leur, Netherlands
+31 76 508 5300
www.parker.com/dhfn

Hiross Zander

Padova Business Unit
Padova, Italy
+39 049 9712 111
www.parker.com/hzd

Hiross Zander

Essen Business Unit
Essen, Germany
+49 2054 9340
www.parker.com/hzd

Engine Filtration & Water Purification

Racor

Dewsbury, England
+44 (0) 1924 487 000
www.parker.com/rfde

Racor Research & Development

Stuttgart, Germany
+49 (0)711 7071 290-10
www.parker.com/rfde

Hydraulic Filtration

Hydraulic Filter

Arnhem, Holland
+31 26 3760376
www.parker.com/hfde

Urjala Operation

Urjala, Finland
+358 20 753 2500
www.parker.com/hfde

Condition Monitoring Center

Norfolk, England
+44 1842 763 299
www.parker.com/hfde

Process Filtration

domnick hunter Process Filtration

Birtley, England
+44 (0) 191 410 5121
www.parker.com/processfiltration

Asia Pacific

Australia

Castle Hill, Australia
+61 2 9634 7777
www.parker.com/australia

China

Shanghai, China
+86 21 5031 2525
www.parker.com/china

India

Navi Mumbai, India
+91 22 651 370 8185
www.parker.com/india

Japan

Tokyo, Japan
+81 45 870 1522
www.parker.com/japan

Korea

Hwaseon-City
+82 31 359 0852
www.parker.com/korea

Singapore

Jurong Town, Singapore
+65 6887 6300
www.parker.com/singapore

Thailand

Bangkok, Thailand
+66 2 186 7000
www.parker.com/thailand

Latin America

Parker Comercio Ltda.

Filtration Division

Sao Paulo, Brazil
+55 12 4009 3500
www.parker.com/br

Pan American Division

Miami, FL
305 470 8800
www.parker.com/panam

Africa

Aeroport Kempton Park, South Africa

+27 11 9610700
www.parker.com/africa

© 2010, 2012 Parker Hannifin Corporation. Product names are trademarks or registered trademarks of their respective companies.

Bulletin AGSLCMS 9/12



Parker Hannifin Corporation
Filtration and Separation Division

242 Neck Road
Haverhill, MA 01835
T 978 858 0505, F 978 556 7501
1 800 343 4048 (US & Canada)
www.labgasgenerators.com



ISO 9001
ISO 14001